




2025 ENT CATALOG

www.shark-instruments.com




ENT

Middle Ear Instruments	01-02
Picks and Scoops	03
Hooks and Probes	04
Retractors	05
Suction Tubes	06
Speculum	07-08
Aural Forceps	09-10
Artery Forceps	11
Tonsil, Turbinate and Septum Forceps	12
Dressing Forceps	13
Other Forceps	14
Mouth Gags	15
Tongue Plates and Depressors	16
Punches	17-20
Rongeurs	21-22
Scissors	23-24
Gouges, Knives and Elevators	25
Curettes and Myringotomes	26
Elevators and Dissectors	27
Miscellaneous	28-29

MIDDLE EAR INSTRUMENTS

Ironside Instruments


160mm long, black finish



Code	Description
SI-ENT/200a	Knife 1, vertical incision
SI-ENT/200b	Knife 2, horizontal incision
SI-ENT/200c	Elevator
SI-ENT/200d	Needle point
SI-ENT/200e	Hook, superior elevation, angle 45°
SI-ENT/200f	Hook, inferior elevation, angle 45°
SI-ENT/200g	Hook, anterior elevation, angle 90° right
SI-ENT/200h	Hook, anterior elevation, angle 90° left
SI-ENT/200i	Measuring device

Rosen Instruments


160mm long, black finish



Code	Description
SI-ENT/202a	Curette, tympanic ring, small
SI-ENT/202b	Curette, tympanic ring, large
SI-ENT/202c	Probe, footplate movement test, fine
SI-ENT/202d	Mobiliser, slight curve
SI-ENT/202e	Mobiliser, full curve
SI-ENT/202f	Scoop, size 000
SI-ENT/202g	Scoop, size 00
SI-ENT/202h	Scoop, size 0
SI-ENT/202i	Scoop, size 1
SI-ENT/202j	Fenestrator, oval window

Schuknecht Instruments


160mm long, black finish



Code	Description
SI-ENT/201a	Hook upward, 0.4mm
SI-ENT/201b	Hook downward, 0.4mm
SI-ENT/201c	Hook upward, 0.7mm
SI-ENT/201d	Hook downward, 0.7mm
SI-ENT/201e	Elevator, oval
SI-ENT/201f	Hook, right-angled upward
SI-ENT/201g	Pick, straight
SI-ENT/201h	Spatula

Beale Instruments

165mm long, satin finish




Code	Description
SI-ENT/204a	Tympanoplasty probe, fine blunt end
SI-ENT/204b	Raspatory, slight curve
SI-ENT/204c	Raspatory, slight curve, small
SI-ENT/204d	Rugine, double ended, sharp
SI-ENT/204e	Rugine, double ended
SI-ENT/204f	Spatula, graft smoothing

MIDDLE EAR INSTRUMENTS

Beale Instruments


170mm long, satin finish



Code	Description
SI-ENT/203a	Knife, Stapedius tendon cutting
SI-ENT/203b	Raspatory, mucous removal
SI-ENT/203c	Osteotome, footplate
SI-ENT/203d	Perforator, elevation fixed footplate
SI-ENT/203e	Hook, removal of footplate, angled downward
SI-ENT/203f	Hook, removal of footplate, angled upward
SI-ENT/203g	Hook, bone removal from oval window, angled left
SI-ENT/203h	Applicator, vein graft for right ear
SI-ENT/203i	Applicator, vein graft for left ear
SI-ENT/203j	Curette, removal bony rim
SI-ENT/203k	Curette

Hough Instruments


165mm long, stainless steel



Code	Description
SI-ENT/205a	Excavator, "whirlybird", right, large
SI-ENT/205b	Excavator, "whirlybird", left, large
SI-ENT/205c	Excavator, "whirlybird", right, small
SI-ENT/205d	Excavator, "whirlybird", left, small
SI-ENT/205e	Pick, stapedectomy footplate, angled 90°, 0.33mm
SI-ENT/205f	Pick, stapedectomy footplate, angled 45°, 0.5mm
SI-ENT/205g	Pick, stapedectomy footplate, angled 15°, 0.5mm
SI-ENT/205h	Hook, "hoe", excavator, angled 90°
SI-ENT/205i	Hook, "hoe", excavator, angled 45°
SI-ENT/205j	Pick, curved

Cawthorne Hooks

165mm long, black finish



Code	Description
SI-ENT/206a	Slight curve
SI-ENT/206b	Medium curve
SI-ENT/206c	Full curve
SI-ENT/206d	Double curve
SI-ENT/206e	90° angle, 0.5mm
SI-ENT/206f	90° angle, 1.0mm
SI-ENT/206g	45° angle, 0.5mm
SI-ENT/206h	Retrograde
SI-ENT/206i	Straight

PICKS AND SCOOPS

Musgrove Cerumen Pick and Wool Carrier


Double ended



Code	Description
SI-ENT/024	178mm long

Formby Cerumen Scoop and Hook

Double ended



Code	Description
SI-ENT/026	178mm long

Hovell Loop


Fine, for removing hard wax



Code	Description
SI-ENT/027	152mm long

Lempert Scoop

Oval, 200mm long



Code	Description
SI-ENT/066a	2mm x 2.5mm
SI-ENT/066b	2mm x 3.5mm
SI-ENT/066c	2mm x 4mm
SI-ENT/066d	2mm x 4.5mm
SI-ENT/066e	3mm x 5.5mm
SI-ENT/066f	3mm x 6mm

House Curette

Double ended, 180mm long



Code	Description
SI-SH/5213	Straight
SI-SH/5213a	Angled

Balance Scoop and Curette



Code	Description
SI-ENT/038	230mm long

HOOKS AND PROBES

Hartmann Hook

Sharp, straight shaft



Code	Description
SI-ENT/021	140mm long


Gardiner Brown Hook and Imray Loop

Double ended



Code	Description
SI-ENT/025	127mm long

Cawthorne Cerumen Hook



Code	Description
SI-ENT/023	143mm long

St Bart's Hook

Blunt, straight handle



Code	Description
SI-ENT/022	178mm long

Tracheal Hook

Sharp



Code	Description
SI-ENT/SH.5672	165mm long

McIndoe Skin Hook

Fine, with flat handle



Code	Description
SI-ENT/057	191mm long

Tracheal Single Hook Retractor

Sharp



Code	Description
SI-GI/SH5671	152mm long

Cawthorne Probe

Silver, flattened end



Code	Description
SI-ENT/019	140mm long

Tracheal Double Hook

Blunt



Code	Description
SI-ENT/SH.5673	105mm long

Jobson Horne Probe

Double ended, serrated tip and smooth ring




Code	Description
SI-ENT/020a	140mm long
SI-ENT/020b	178mm long

RETRACTORS


Mollison Retractor

Self retaining, curved



Code	Description
SI-ENT/034a	4-4 teeth, adult, 165mm long
SI-ENT/034b	2-2 teeth, child, 133mm long


Wallace Modified Plester Retractor



Code	Description
SI-ENT/062	2 x 2 teeth, 140mm long

Plester-Jansen Retractor


Solid blade one side,
140mm long



Code	Description
SI-Ent/061A	Left
SI-Ent/061B	Right
SI-Ent/061C	3 x 3 Teeth

Shah Intra Aural Retractor

Satin finish
152mm long



Code	Description
SI-ENT/060a	4-4 teeth, adult, 165mm long
SI-ENT/060b	Third blade
SI-ENT/060c	Intra aural blade
SI-ENT/060d	Post aural blade

Hurd Pillar Retractor and Tonsil Dissector



Code	Description
SI-ENT/147	197mm long

Mollison Anterior Pillar Retractor and Tonsil Enucleator


Semi sharp



Code	Description
SI-ENT/148	210mm long

SUCTION TUBES

Ipswich Pattern Sucker



Code	Description
SI-ENT/050a	267mm long

Lempert Suction and Irrigation Tubes


With tucker valve and stilette for cleaning stainless steel, 203mm



Code	Description
SI-ENT/044	6F.G. 2mm
SI-ENT/044a	9F.G. 3mm

Formby Suction Tube

Angled, conical, tip size 10 F.G. stainless steel



Code	Description
SI-ENT/042	210mm long

Magill Suction Tubes


Set of four with conical ends, one curved, three angled, handle, stainless steel, 210mm



Code	Description
SI-ENT/045	Various sizes
SI-ENT/046a	3.0mm ø
SI-ENT/046b	3.5mm ø
SI-ENT/046c	4.0mm ø

Barron Suction Tube


Small size, 4 F.G.
Not for use with Zoellner ends



Code	Description
SI-ENT/041	229mm long

Zoellner Suction Tube

With stilette, stainless steel




Code	Description
SI-ENT/047	170mm long, 2mm ø




House Suction and Irrigation Tubes

Stainless steel, 152mm long



Code	Description
SI-ENT/043	8 F.G. x 6 F.G.
SI-ENT/043a	10 F.G. x 7 F.G.
SI-ENT/043b	12 F.G. x 8 F.G.

Rowbotham Suction Tube



Code	Description
SI-ENT/169	286mm long

SPECULUM

Rosen Slotted Speculum

Set of four, round ends, black finish
38mm long

Code	Description
SI-ENT/049j	Set of 4
SI-ENT/049f	4mm
SI-ENT/049g	5mm
SI-ENT/049h	6mm
SI-ENT/049i	7mm



Rosen Aural Speculum

Round, 40mm long

Code	Description
SI-ENT/049a	3mm
SI-ENT/049b	4mm
SI-ENT/049c	5mm
SI-ENT/049d	6mm
SI-ENT/049e	7mm



Shea Speculum

Set of five, oblique ends, black finish, 38mm long


Code	Description
SI-ENT/050	Set of 5
SI-ENT/051	3mm
SI-ENT/051a	4mm
SI-ENT/051b	5mm
SI-ENT/051c	6mm
SI-ENT/051d	7mm



Tumarkin Speculum

Set of eight, tapered slot, black finish
50mm long

Code	Description
SI-ENT/052	Set of 8
SI-ENT/053a	4.5mm ø
SI-ENT/053b	6.5mm ø
SI-ENT/053c	7.0mm ø
SI-ENT/053d	7.5mm ø
SI-ENT/053e	8.0mm ø
SI-ENT/053f	8.5mm ø
SI-ENT/053g	9.0mm ø
SI-ENT/053h	11.0mm ø




Heath Speculum

Set of four, black finish
50mm long

Code	Description
SI-ENT/207e	Set of 4
SI-ENT/207a	4.0mm ø
SI-ENT/207b	5.0mm ø
SI-ENT/207c	6.0mm ø
SI-ENT/207d	7.0mm ø



Peters Seigle Speculum




Code	Description
SI-ENT/366	Peters ear speculum black chrome finish set of 5, sizes 0 to 4
SI-ENT/189i	2.2mm x 3.2mm
SI-ENT/189j	2.8mm x 4.2mm
SI-ENT/189k	4.3mm x 5.8mm
SI-ENT/189l	5.2mm x 7.6mm
SI-ENT/189m	6.0mm x 7.6mm
SI-ENT/189d	6mm
SI-ENT/189e	7mm
SI-ENT/189f	8mm
SI-ENT/189g	9mm
SI-ENT/189h	10mm

SPECULUM

St Clair Thompson Nasal Speculum

Code	Description
SI-ENT/072	38mm long
SI-ENT/072a	51mm long
SI-ENT/072b	64mm long
SI-ENT/072c	89mm long



Thudichum Nasal Speculum

Set of seven, sizes 1-7,
70mm long

Code	Description
SI-ENT/073a	Standard
SI-ENT/073b	Insulated



Standard

Code	Description
SI-ENT/073a	Size 1 - 12mm blade
SI-ENT/073a	Size 2 - 14mm blade
SI-ENT/073a	Size 3 - 16mm blade
SI-ENT/073a	Size 4 - 18mm blade
SI-ENT/073a	Size 5 - 20mm blade
SI-ENT/073a	Size 6 - 22mm blade
SI-ENT/073a	Size 7 - 24mm blade
SI-ENT/073a	Size 8 - 27mm blade


Insulated

Code	Description
SI-ENT/073b	Size 1 - 12mm blade
SI-ENT/073b	Size 2 - 14mm blade
SI-ENT/073b	Size 3 - 16mm blade
SI-ENT/073b	Size 4 - 18mm blade
SI-ENT/073b	Size 5 - 20mm blade
SI-ENT/073b	Size 6 - 22mm blade
SI-ENT/073b	Size 7 - 24mm blade
SI-ENT/073b	Size 8 - 27mm blade

Killian Nasal Speculum

With Stop


Code	Description
SI-ENT/081a	51mm blade
SI-ENT/081b	64mm blade
SI-ENT/081c	76mm blade
SI-ENT/081d	89mm blade



Without Stop


Code	Description
SI-ENT/082b	76mm blade

Lempert Endaural Speculum



Code	Description
SI-ENT/060e	133mm long


Pilcher Nasal Speculum



Code	Description
SI-ENT/071	152mm

Cottle Nasal Speculum

127mm long




Code	Description
SI-ENT/192c	20mm blade
SI-ENT/192b	32mm blade
SI-ENT/192	41mm blade
SI-ENT/192a	60mm blade


AURAL FORCEPS



Triangular Jaw




Smooth Jaw



Serrated Jaw




Tenaculum Jaw



Cup Jaw



Grunwald Punch



Henckel Jaw

Beale Aural Forceps

Micro, black finish, triangular flat jaws, serrated



Code	Description
SI-ENT/183	70mm to shoulder

Hough Aural Forceps

Micro, extra fine oval cup jaws, black finish



Code	Description
SI-ENT/059	70mm to shoulder

Cawthorne Aural Forceps

Micro, black finish, extra fine
70mm to shoulder



Code	Description
SI-ENT/182	Smooth jaws
SI-ENT/199	Serrated jaws

Ormerod Aural Forceps


Micro, extra fine, black finish
70mm to shoulder



Code	Description
SI-ENT/184a	Serrated, angled upwards
SI-ENT/184b	1mm cup jaw
SI-ENT/184c	Elongated cup jaw, 0.5 x 4mm

Hartmann Aural Forceps

Micro, black finish
75mm to shoulder



Code	Description
SI-ENT/058	Serrated jaws
SI-ENT/058a	Tenaculum jaws
SI-ENT/058b	Cup shape jaws
SI-ENT/058c	Curved left, serrated
SI-ENT/058d	Curved right, serrated
SI-ENT/058e	Angled up, serrated
SI-ENT/058f	Serrated jaws, satin finish
SI-ENT/058s	Pre-sterile - single use only

Shea Aural Forceps

Micro, oval cup jaw, black finish
70mm to shoulder




Code	Description
SI-ENT/069	Standard
SI-ENT/069a	Curved to left
SI-ENT/069b	Curved to right
SI-ENT/069c	Curved upward

AURAL FORCEPS

Hartmann Aural Forceps


Standard, satin finish
75mm to shoulder



Code	Description
SI-ENT/181	Serrated jaws
SI-ENT/181a	Tenaculum jaws
SI-ENT/181b	Cup shape jaws
SI-ENT/181c	Grunwald punch action

Heath Aural Forceps


Light, with lifting discs, serrated jaws



Code	Description
SI-ENT/175	83mm to shoulder

Henckel Aural Forceps


Fine fenestrated scoop blades



Code	Description
SI-ENT/176	75mm to shoulder

Struycken Cutting Forceps

75mm to shoulder




Code	Description
SI-ENT/178	Fenestrated under jaw

ARTERY FORCEPS

Birkett Tonsil Artery Forceps


Fine, 178mm long



Code	Description
SI-ENT/151a	Straight
SI-ENT/151b	Curved


Negus Tonsil Artery Forceps

180mm long



Code	Description
SI-ENT/197	Large curve
SI-ENT/197a	Small curve


Denis Browne Tonsil Artery Forceps



Code	Description
SI-ENT/155	216mm long

Sawtell Tonsil Artery Forceps


178mm long



Code	Description
SI-ENT/152	Large curve
SI-ENT/153	Small curve

Fickling Artery Forceps


190mm long



Code	Description
SI-PL/044b	1-2 teeth
SI-PL/044a	Serrated

Scott Tonsil Artery Forceps

Curved to side, screw joint




Code	Description
SI-ENT/154	180mm long

TONSIL, TURBINATE AND SEPTUM FORCEPS

Denis Browne Tonsil Holding Forceps


Angled shanks, no rack



Code	Description
SI-ENT/144	200mm long

Jansen-Middleton Nasal Septum Forceps


Compound action
178mm long



Code	Description
SI-PL/032	5mm jaw


Irwin Moore Nasal Turbinate Forceps

191mm long




Code	Description
SI-ENT/077d	5mm jaw

Walsham Septum Forceps



Code	Description
SI-PL/030a	Straight, 215mm long
SI-PL/030b	Right, 215mm long
SI-PL/030c	Left, 215mm long

Luc Nasal Turbinate Forceps



Heart shaped, 178mm long


Code	Description
SI-ENT/078a	Small 6mm ø
SI-ENT/078b	Medium 7.5mm ø
SI-ENT/078c	Large 9.0mm ø

Oval, 191mm long

Code	Description
SI-ENT/078d	Small 6mm ø
SI-ENT/078e	Medium 7.5mm ø
SI-ENT/078f	Large 9.0mm ø

Hartmann Conchotomes

Circular cutting
102mm to shoulder




Code	Description
SI-ENT/080a	5mm ø
SI-ENT/080b	7mm ø
SI-ENT/080c	9mm ø
SI-ENT/080d	11mm ø

DRESSING FORCEPS

Hartmann Dressing Forceps


Serrated jaw, screw joint



Code	Description
SI-ENT/076	108mm to shoulder

Wilde Dressing Forceps


Angular, serrated jaws



Code	Description
SI-ENT/075	102mm to shoulder

Tilley Dressing Forceps


Straight, fine, serrated ends



Code	Description
SI-ENT/077	2mm wide, 76mm to shoulder

Irwin Moore Thyroid Clamp


Thin jaws, angled to side



Code	Description
SI-ENT/170	165mm long

OTHER FORCEPS


Trousseau Tracheal Dilator



Code	Description
SI-ENT/172	Adult 140mm long
SI-ENT/172a	Child 125mm long

Wilson Bone Nibbling Forceps


Fine, curved, satin finish
152mm long



Code	Description
SI-ENT/054	Downward cutting - 3mm jaws
SI-ENT/055	Upward cutting - 3mm jaws


Cawthorne Nibbling Forceps

Narrow, 3mm jaws, wide handle



Code	Description
SI-ENT/068	159mm long


Denis Browne Collar Crushing Forceps



Code	Description
SI-4062	140mm long

Lempert Nibbling Forceps


152mm long



Code	Description
SI-ENT/067a	Straight, 2.5mm jaws
SI-ENT/067b	Angled, 3.0mm jaws

MOUTH GAGS


Doyen Mouth Gag



Code	Description
SI-ENT/88d	Adult
SI-ENT/88e	Child

Boyle-Davis Mouth Gag


Frame only



Code	Description
SI-ENT/083a	Adult
SI-ENT/083b	Child


Kilner Dott Mouth Gag

Frame only




Code	Description
SI-ENT/087a	Adult
SI-ENT/088e	Set, including frame, small, medium, and large blades

Fergusson Ackland Mouth Gag




Code	Description
SI-ENT/209	Adult, 210mm
SI-ENT/209a	Child, 155mm

Mason Mouth Gag



Code	Description
SI-ENT/210	Adult, 210mm

Draffin Bipod



Adult

Code	Description
SI-ENT/083	Adult, 4 hole x 480mm long half circle end for use with a drape
SI-ENT/083E	Adult, 4 hole x 480mm long with round ends to be used with Magauran plate
SI-ENT/083E	Adult, 4 hole x 480mm long with round ends to be used with Magauran plate

Child


Code	Description
SI-ENT/083C	Child, 3 hole x 425mm long half circle end for use with a drape
SI-ENT/083D	Child, 3 hole x 425mm long with round ends to be used with Magauran plate
SI-ENT/083F	Child, 4 hole x 425mm long with round ends to be used with Magauran plate Alder Hey Children's Hospital pattern
SI-ENT/083D	Child, 3 hole x 425mm long with round ends to be used with Magauran plate
SI-ENT/083F	Child, 4 hole x 425mm long with round ends to be used with Magauran plate Alder Hey Children's Hospital pattern
SI-5104	Magauran plate 480mm - for use with Draffin bipods

Magauran Plate

Code	Description
SI-5104	480mm - for use with Draffin bipods

TONGUE PLATES AND DEPRESSORS

Doughty Tongue Plate



Plain

Code	Description
SI-ENT/084a	65mm - 1.6mm-1.9mm slot
SI-ENT/084b	75mm - 1.6mm-1.9mm slot
SI-ENT/084c	90mm - 4mm-4.5mm slot
SI-ENT/084d	100mm - 4mm-4.5mm slot
SI-ENT/084e	115mm - 4mm-4.5mm slot

Slotted


Code	Description
SI-ENT/085a	65mm - 1.6mm-1.9mm slot
SI-ENT/085b	75mm - 1.6mm-1.9mm slot
SI-ENT/085c	90mm - 4mm-4.5mm slot
SI-ENT/085d	100mm - 4mm-4.5mm slot
SI-ENT/085e	115mm - 4mm-4.5mm slot

With anaesthetic tube

Code	Description
SI-ENT/086a	65mm - 1.6mm-1.9mm slot
SI-ENT/086b	75mm - 1.6mm-1.9mm slot
SI-ENT/086c	90mm - 4mm-4.5mm slot
SI-ENT/086d	100mm - 4mm-4.5mm slot
SI-ENT/086e	115mm - 4mm-4.5mm slot

Kilner Dott Tongue Plate

With endotracheal tube groove



Code	Description
SI-ENT/088a	Small
SI-ENT/088b	Medium
SI-ENT/088c	Large
SI-ENT/088d	Extra small


Tongue Depressor

Flat stainless steel



Code	Description
SI-ENT/143	184mm long


Lack Tongue Depressor



Code	Description
SI-ENT/142a	Child, 13mm wide
SI-ENT/142b	Adolescent, 19mm wide
SI-ENT/142c	Adult, 22mm wide

PUNCHES


Nasal Cutting Forceps



Myles

Through cutting, 120mm to shoulder


Code	Description
SI-ENT/SH5889	12.5 x 3.5mm jaw



Gruenwald-Henke

Through cutting, 120mm to shoulder


Code	Description
SI-ENT/SH5281	7 x 2mm jaw
SI-ENT/SH5282	9 x 2.5mm jaw



Struycken


Through cutting, 120mm to shoulder

Code	Description
SI-ENT/SH5280	11 x 1mm jaw



Gruenwald Nasal Punch Forceps


108mm to shoulder



Code	Description
SI-ENT/077a	A fenestrated jaw, 3.5mm
SI-ENT/077b	Fenestrated lower jaw, narrow, 2.5mm
SI-ENT/077c	Fenestrated and tapered lower jaw, 2-3mm

Seymour Jones Antrum Punch Forceps


159mm long



Code	Description
SI-ENT/132	4mm jaw

Backbiter


180mm long



Code	Description
SI/4765b	7 x 3mm jaw

Backbiter

200mm long




Code	Description
SI-ENT/5392	7 x 2mm jaw

PUNCHES

Stammberger Rhinoforcep Antrium Punch


7mm x 3mm, 130mm long



Code	Description
SI-ENT/234	Left back cut
SI-ENT/234a	Right back cut

Tilley Henckel Ethmoid Punch


Graduated shank,
114mm to shoulder



Code	Description
SI-ENT/138a	4.5mm
SI-ENT/138b	5.5mm
SI-ENT/138c	6.5mm
SI-ENT/138d	8mm

Ostrom Antrum Punch


Forcep, 178mm to shoulder



Code	Description
SI-ENT/135	3mm jaw

Mackay Modified Stammberger Punch

200mm long



Code	Description
SI-ENT/5392	7 x 2mm jaw

Paterson Biopsy Forceps

Elongated cup, 3mm jaw




Code	Description
SI-ENT/174a	305mm long
SI-ENT/174b	355mm long
SI-ENT/174c	406mm long
SI-ENT/174d	508mm long


PUNCHES

Chevalier Jackson Endoscope Forcep

With crocodile action and serrated jaw



Code	Description
SI-ENT/241	Effective shaft length 300mm
SI-ENT/242	Effective shaft length 400mm
SI-ENT/243	Effective shaft length 500mm
SI-ENT/244	Effective shaft length 600mm




Chevalier Jackson Vocal Nodule Endoscope Forcep

With crocodile action and curved cup



Code	Description
SI-ENT/245	Effective shaft length 300mm
SI-ENT/246	Effective shaft length 400mm
SI-ENT/247	Effective shaft length 500mm
SI-ENT/248	Effective shaft length 600mm




Chevalier Jackson Vocal Nodule Endoscope Forcep

With crocodile action and straight cup




Code	Description
SI-ENT/249	Effective shaft length 300mm
SI-ENT/250	Effective shaft length 400mm
SI-ENT/251	Effective shaft length 500mm
SI-ENT/252	Effective shaft length 600mm



St Bartholomew's Ethmoid Punch

Scoop blades, 3mm bite, 114mm to shoulder




Code	Description
SI-ENT/140a	Straight
SI-ENT/140b	Upward cutting
SI-ENT/140c	Downward cutting

PUNCHES

Ferris Smith Kerrison Punch


Heavy, 5mm bite,
229mm long



Code	Description
SI-ENT/137a	Upward cutting
SI-ENT/137B	Downward cutting

Citelli Antrum Punch


76mm to shoulder



Code	Description
SI-ENT/133a	3mm jaw
SI-ENT/133b	5mm jaw

Hajek Antrum Punch


197mm long



Code	Description
SI-ENT/134	5mm jaw

Hajek Sphenoidal Punch

146mm to shoulder



Upward cutting

Code	Description
SI-ENT/141a	4mm bite
SI-ENT/141b	5mm bite


Downward cutting

Code	Description
SI-ENT/141c	4mm bite
SI-ENT/141d	5mm bite

RONGEURS

Beck Rongeur


Straight, for intervertebral disc



Code	Description
SI-NEU/125	2mm bite, 130mm shaft

Pennybacker Rongeur


120mm shaft



Code	Description
SI-NEU/133	3mm bite
SI-NEU/107	4mm bite

Cushing Rongeur

130mm shaft




Code	Description
SI-NEU/127-2	2mm bite
SI-NEU/127-3	3mm bite
SI-NEU/127-4	4mm bite
SI-NEU/127-5	5mm bite
SI-NEU/127-6	6mm bite



St Bart's Rongeur


114mm shaft



Code	Description
SI-ENT/140a	3mm bite

Cushing Rongeur

130mm shaft



150° upward cutting


Code	Description
SI-NEU/128c	2mm bite
SI-NEU/128	3mm bite
SI-NEU/128a	4mm bite
SI-NEU/128b	5mm bite

150° downward cutting

Code	Description
SI-NEU/130c	2mm bite
SI-NEU/130	3mm bite
SI-NEU/130a	4mm bite
SI-NEU/130b	5mm bite

St Bart's Rongeur

140°, 3mm, 114mm shaft




Code	Description
SI-ENT/140b	Upward cutting
SI-ENT/140ac	Downward cutting

RONGEURS

Mackay Rongeur


140°, straight, 114mm shaft



Code	Description
SI-SH/4249	5.5mm

Weil-Blakesley Rongeur

110mm shaft



140° upward cutting


Code	Description
SI-ENT/194b	3.6mm bite
SI-ENT/194e	4.2mm bite
SI-ENT/194f	4.8mm bite

140° downward cutting

Code	Description
SI-ENT/194g	3.6mm bite
SI-ENT/194h	4.2mm bite
SI-ENT/194i	4.8mm bite

Mackay Rongeur


With suction tube



Code	Description
SI-SH/4249a	3.5mm bite, 114mm shaft

Weil-Blakesley Rongeur


120mm shaft



Code	Description
SI-ENT/193b	3.0mm bite
SI-ENT/193	3.6mm bite
SI-ENT/194	4.2mm bite
SI-ENT/194a	4.8mm bite
SI-ENT/194c	5.8mm bite
SI-ENT/194d	6.6mm bite

Mackay Rongeur

With suction tube, 140° upward cutting




Code	Description
SI-SH/4249b	3.5mm bite, 114mm shaft

SCISSORS

Beckmann Turbinectomy Scissors


Curved



Code	Description
SI-ENT/097	165mm long


Fomon Scissors

Straight, serrated blades, angled handles



Code	Description
SI-ENT/200	145mm long


Scott Stevenson Tonsil Scissors



Code	Description
SI-ENT/159	216mm long

Cottle Scissors


Straight blades, angled handles, blunt/blunt



Code	Description
SI-ENT/OK201	165mm long

Wood Tonsil Scissors


Sharpened outer edges



Code	Description
SI-ENT/156	178mm long

Cottle Bulldog Scissors


Blunt/blunt, 115mm long



Code	Description
SI-ENT/202	115mm long

Wheeler Tonsil Scissors


Triangular form blades, serrated shanks



Code	Description
SI-ENT/157	184mm long

Joseph Scissors


Sharp/sharp, 150mm long



Code	Description
SI-ENT/203	Straight
SI-ENT/204	Curved

Tilley Tonsil Scissors


Curved



Code	Description
SI-ENT/158	190mm long

Heymann Turbinectomy Scissors


Angular



Code	Description
SI-ENT/098	197mm long

SCISSORS


Wilson Tonsil Scissors



Code	Description
SI-ENT/160	203mm long


Cawthorne Scissors

Micro, extra fine, 70mm to shoulder



Code	Description
SI-ENT/186c	Left
SI-ENT/186d	Right
SI-ENT/186e	Straight


Thyroid Scissors



Code	Description
SI-ENT/171	152mm long

Ormorod Aural Scissors


Micro, extra fine, curved upwards



Code	Description
SI-ENT/185	70mm to shoulder

Knapp Scissors


Slight curve, blunt points, 114mm long



Code	Description
SI-ENT/187a	Satin finish
SI-ENT/187b	Black finish

Portmann Aural Scissors

Micro, satin finish 70mm to shoulder




Code	Description
SI-ENT/186a	Curved to right
SI-ENT/186b	Curved to left

GOUGES, KNIVES AND ELEVATORS

Ballenger Nasal Cartilage Knife

Revolving blade



Straight, 203mm long

Code	Description
SI-ENT/096a	6mm wide
SI-ENT/096b	8mm wide

Cranked, 197mm long

Code	Description
SI-ENT/096c	6mm wide

Sickle Knife



Code	Description
SI-BUPNS002	188mm long
SI-ENT/271	160mm long

Jenkins Gouge

140mm long



Code	Description
SI-ENT/037	1mm
SI-ENT/037a	2mm
SI-ENT/037b	3mm
SI-ENT/037c	4mm
SI-ENT/037d	5mm
SI-ENT/037e	6mm
SI-ENT/037f	7mm

Code	Description
SI-ENT/037g	8mm
SI-ENT/037h	9mm
SI-ENT/037i	10mm
SI-ENT/037j	11mm
SI-ENT/037k	12mm
SI-ENT/037l	13mm
SI-ENT/037m	14mm

Killian Nasal Gouge

127mm long



Code	Description
SI-ENT/093	4mm

Guy's Hospital Nasal Gouge

Shallow groove and flattened back, 127mm long



Code	Description
SI-ENT/094	4mm

Tilley Nasal Gouge

V shaped, 127mm long



Code	Description
SI-ENT/095	4mm

Cottle Elevator

Double ended, black chrome finish
200mm long



Code
SI-ENT/224

Hill Elevator

Double ended, curved to side, rugine end



Code	Description
SI-ENT/091	5mm wide, 222mm long

Tumarkin Atticotomy Elevator

Double ended




Code	Description
SI-ENT/036	191mm long

CURETTES AND MYRINGOTOMES

Buck Ear Curette

Angled, blunt, 146mm long



Code	Description
SI-ENT/028	2mm
SI-ENT/028a	2.5mm
SI-ENT/028b	3mm
SI-ENT/028c	3.5mm
SI-ENT/028d	4mm

Nasal Curette


Straight, 185mm long



Code	Description
SI-ENT/628702	5 x 2mm

Beckmann Adenoid Curette


216mm long



Code	Description
SI-ENT/167a	8mm blade
SI-ENT/167b	10mm blade
SI-ENT/167c	12mm blade
SI-ENT/167d	14mm blade
SI-ENT/167e	16mm blade
SI-ENT/167f	18mm blade

Agnew (Politzer) Myringotome

190mm long



Code	Description
SI-ENT/031	Shaft length 100mm


Trautmann Myringotome

178mm long



St Clair Thomson Adenoid Curette

With cage, 235mm long




Code	Description
SI-ENT/168a	8mm
SI-ENT/168b	10mm
SI-ENT/168c	12mm
SI-ENT/168d	14mm
SI-ENT/168e	16mm
SI-ENT/168f	18mm

Code	Description
SI-ENT/032	Upward cutting
SI-ENT/033	Downward cutting

ELEVATORS AND DISSECTORS

Freer Instruments

Single ended, 190mm long



Knife

Code	Description
SI-ENT/065a	Curved to side
SI-ENT/065c	Rounded edge, 8mm blade
SI-ENT/065d	Flat rounded blade
SI-ENT/065e	Rounded blade, curved on flat

Elevator

Code	Description
SI-ENT/065f	Dull edged
SI-ENT/065h	Sharp edged

Chisel

Code
SI-ENT/065i

St Clair Thompson Dissector

159mm



Code	Description
SI-ENT/089a	Sharp
SI-ENT/089b	Blunt

Negus Aspiring Dissector


With 6mm olive, 210mm long



Code	Description
SI-ENT/090a	-
SI-ENT/090b	Insulated

Freer Separator and Elevator

Double ended



Code	Description
SI-ENT/065b	Blunt/sharp, 180mm long
SI-ENT/065g	Sharp, knife curved to side, 180mm long
SI-ENT/065j	Blunt/sharp, 200mm long
SI-ENT/065j-sal	Blunt/sharp, 188mm long

Syme Tonsil Dissector

Double ended, one end blunt, one semi sharp



Code	Description
SI-ENT/149	170mm long

Gwynne-Evans Tonsil Dissector

Double ended, with serrated end and blunt end
197mm long




Code	Description
SI-ENT/146a	6mm wide
SI-ENT/146b	8mm wide

MISCELLANEOUS

Woakes Eustachian Catheters

Silver, Luer, 140mm



Code	Description
SI-ENT/015	7 F.G.
SI-ENT/015a	8 F.G.
SI-ENT/015b	9 F.G.

Code	Description
SI-ENT/015c	10 F.G.
SI-ENT/015d	11 F.G.

Kramer Eustachian Catheters

Silver, Luer, 140mm




Code	Description
SI-ENT/014	7 F.G.
SI-ENT/014a	8 F.G.
SI-ENT/014b	9 F.G.

Code	Description
SI-ENT/014c	10 F.G.
SI-ENT/014d	11 F.G.

Myles Naso-Antral Perforator


171mm long



Code	Description
SI-ENT/103a	16 F.G.
SI-ENT/103b	19 F.G.

Tilley Antrum Burrs


171mm long



Code	Description
SI-ENT/131a	6mm ø
SI-ENT/131b	7mm ø
SI-ENT/131c	8mm ø

Laryngeal Mirror

Mirror only, 130mm long




Code	Description
SI-ENT/191-6	6mm ø
SI-ENT/191-8	8mm ø
SI-ENT/191-10	10mm ø
SI-ENT/191-12	12mm ø
SI-ENT/191-14	14mm ø
SI-ENT/191-16	16mm ø
SI-ENT/191-18	18mm ø

Code	Description
SI-ENT/191-20	20mm ø
SI-ENT/191-22	22mm ø
SI-ENT/191-24	24mm ø
SI-ENT/191-26	26mm ø
SI-ENT/191-28	28mm ø
SI-ENT/191-30	30mm ø

Laryngeal Mirror

Complete with metal handle, 260mm long



Code	Description
SI-ENT/192-6	6mm ø
SI-ENT/192-8	8mm ø
SI-ENT/192-10	10mm ø
SI-ENT/192-12	12mm ø
SI-ENT/192-14	14mm ø
SI-ENT/192-16	16mm ø
SI-ENT/192-18	18mm ø

Code	Description
SI-ENT/192-20	20mm ø
SI-ENT/192-22	22mm ø
SI-ENT/192-24	24mm ø
SI-ENT/192-26	26mm ø
SI-ENT/192-28	28mm ø
SI-ENT/192-30	30mm ø

Mitchell Trimmer



Code	Description
SI-ENT/013	165mm

MISCELLANEOUS

Lake Seeker

Double ended, for vestibule



Code	Description
SI-ENT/040	184mm long


Negus Pusher

"Knot tier", metal handle




Code	Description
SI-ENT/165	203mm long

Lempert Rugine



Code	Description
SI-ENT/064a	Narrow, 156mm long
SI-ENT/065a	Broad, 171mm long

Stanford Cade Insufflator



Code	Description
SI-ENT/188	175mm long

Irwin-Moore Tonsil Needle

203mm long



Code	Description
SI-ENT/164a	Curved to right
SI-ENT/164a	Curved to left

Tubing Mount

For use with catheters, Luer



Code	Description
SI-ENT/016	25mm



Contact Information



www.shark-instruments.com



info@shark-instruments.com



052-4598813



Mohallah Chah Tarkhana Nargate Street
Jinnah Park, Colony no.2, Sialkot.