




2025 GENERAL CATALOG

www.shark-instruments.com




General Instruments

Artery Forceps	01
Clamps	04
Spring Forceps	05
Tissue Forceps	07
Sponge Holding and Sterilising Forceps	09
Other Forceps	10
Scissors	12
Needle Holders	16
Combined Scissor and Needle Holders	20
Bag and Towel Clips	21
Retractors	22
Dissectors and Probes	30
Miscellaneous	31
Hollow Ware	32
Silicone Matting	33
Sterilising Baskets	34
Optical Fibre Cables	36

ARTERY FORCEPS

Adson Artery Forceps

Box joint, 180mm long



Serrated jaws

Code	Description
SI-GI/034	Straight
SI-GI/034a	Curved




1-2 teeth

Code	Description
SI-GI/034b	Straight
SI-GI/034c	Curved

Cushing Artery Forceps


Box joint, fine, half serrated jaws, 140mm long



Code	Description
SI-GI/001	Straight
SI-GI/002	Curved

Cairns Artery Forceps

Box joint, half serrated jaws
145mm long



Serrated jaws

Code	Description
SI-GI/415	Straight




1-2 teeth

Code	Description
SI-GI/020	Curved
SI-GI/021	Curved to side

Dunhill Artery Forceps

Box joint, fine serrations



Straight


Code	Description
SI-GI/065	135mm long
SI-GI/065a	190mm long

Curved

Code	Description
SI-GI/043	135mm long
SI-GI/468	180mm long

Crile Artery Forceps

Box joint, serrated jaws



Straight

Code	Description
SI-GI/014	140mm long
SI-GI/014a	160mm long




Curved

Code	Description
SI-GI/015	140mm long
SI-GI/015a	160mm long

Lawson Tait Artery Forcep

Box joint, straight jaw




Code	Description
SI-GI/SH5924	140mm long

ARTERY FORCEPS

Halstead Mosquito Forceps

Standard pattern, haemostatic, box joint



Serrated jaw

Code	Description
SI-GI/007	Straight 125mm long
SI-GI/008	Curved 125mm long
SI-GI/009	Fully curved 125mm long

1-2 teeth


Code	Description
SI-GI/010	Straight 125mm long
SI-GI/012	Straight 190mm long
SI-GI/011	Curved 125mm long
SI-GI/013	Curved 190mm long

Micro, 100mm long

Code	Description
SI-GI/003	Straight, serrated
SI-GI/004	Curved, serrated
SI-GI/005	Straight, 1-2 teeth
SI-GI/006	Curved, 1-2 teeth

Heiss Artery Forceps

Box joint, 205mm long



Serrated jaws

Code	Description
SI-GI/016	Straight
SI-GI/017	Curved

1-2 teeth

Code	Description
SI-GI/016a	Straight
SI-GI/017b	Curved

Kelly Artery Forceps

Box joint, half serrated jaws




Straight

Code	Description
SI-GI/022	140mm long
SI-GI/023	160mm long

Curved

Code	Description
SI-GI/024	140mm long
SI-GI/025	160mm long

Fraser Kelly Artery Forceps Curved




Straight

Code	Description
SI-GI/468a	180mm

Kocher Artery Forceps

Box joint, 1-2 teeth



Straight

Code	Description
SI-GI/052	130mm long
SI-GI/053	150mm long
SI-GI/054	180mm long
SI-GI/055	200mm long
SI-GI/056	230mm long


Curved

Code	Description
SI-GI/057	150mm long
SI-GI/058	180mm long
SI-GI/059	200mm long
SI-GI/060	230mm long

ARTERY FORCEPS

Moynihan Artery Forceps

Box joint, half-serrated jaws



Straight


Code	Description
SI-GI/018	145mm long
SI-GI/068	180mm long

Curved

Code	Description
SI-GI/019	145mm long
SI-GI/069	180mm long

Rochester-Pean Artery Forceps

Box joint



Straight


Code	Description
SI-GI/044	165mm long
SI-GI/045	190mm long
SI-GI/046	205mm long
SI-GI/047	240mm long

Curved

Code	Description
SI-GI/048	165mm long
SI-GI/049	180mm long
SI-GI/050	205mm long
SI-GI/051	240mm long

Roberts Artery Forceps


Screw joint, curved on flat, serrated



Code	Description
SI-GI/070	230mm long

Spencer Wells Artery Forceps

Box joint, serrated jaws



Straight

Code	Description
SI-GI/035	127mm long
SI-GI/036	152mm long
SI-GI/037	177mm long
SI-GI/038	203mm long


Curved

Code	Description
SI-GI/039	127mm long
SI-GI/040	152mm long
SI-GI/041	177mm long
SI-GI/042	203mm long

OTHER SIZES AVAILABLE

Dandy Artery Forceps


Box joint, curved in plane



Code	Description
SI-GI/076	140mm long

Carmalt Angiotribe Forceps

Box joint, atraumatic, light pattern



Straight


Code	Description
SI-GI/062	165mm long
SI-GI/063	200mm long
SI-GI/064	230mm long

Curved

Code	Description
SI-GI/061	165mm long

CLAMPS

Dieffenbach Bulldog Clamp



Straight, serrated

Code	Description
SI-GI/112	1.2mm fine jaw, 40mm long
SI-GI/113	1.2mm fine jaw, 45mm long
SI-GI/108	1.8mm jaw, 30mm long
SI-GI/109	1.8mm jaw, 35mm long
SI-GI/110	2.5mm jaw, 40mm long
SI-GI/111	2.5mm jaw, 45mm long
SI-GI/107	6mm heavy jaw, 60mm long
SI-GI/455	20mm jaw, 75mm long
SI-GI/455a	30mm jaw, 85mm long
SI-GI/455b	45mm jaw, 105mm long
SI-GI/455c	65mm jaw, 120mm long

Straight, plain

Code	Description
SI-GI/129	1.2mm fine jaw, 40mm long
SI-GI/131	1.2mm fine jaw, 45mm long
SI-GI/121	1.8mm jaw, 30mm long
SI-GI/123	1.8mm jaw, 35mm long
SI-GI/125	2.5mm jaw, 40mm long
SI-GI/127	2.5mm jaw, 45mm long

Curved, serrated

Code	Description
SI-GI/119	1.2mm fine jaw, 40mm long
SI-GI/120	1.2mm fine jaw, 45mm long
SI-GI/115	1.8mm jaw, 30mm long
SI-GI/116	1.8mm jaw, 35mm long
SI-GI/117	2.5mm jaw, 40mm long
SI-GI/118	2.5mm jaw, 45mm long
SI-GI/114	6mm heavy jaw, 60mm long
SI-GI/456	20mm jaw, 75mm long
SI-GI/456a	30mm jaw, 85mm long
SI-GI/456b	45mm jaw, 105mm long
SI-GI/456c	65mm jaw, 120mm long

Curved, plain

Code	Description
SI-GI/130	1.2mm fine jaw, 40mm long
SI-GI/132	1.2mm fine jaw, 45mm long
SI-GI/122	1.8mm jaw, 30mm long
SI-GI/124	1.8mm jaw, 35mm long
SI-GI/126	2.5mm jaw, 40mm long
SI-GI/128	2.5mm jaw, 45mm long

Blalock Bulldog Clamp



Serrated, with hole in handle

Straight

Code	Description
SI-GI/480	11mm jaw, 38mm long
SI-GI/482	14mm jaw, 48mm long
SI-GI/493	17mm jaw, 60mm long
SI-GI/481	27mm jaw, 70mm long

Curved

Code	Description
SI-GI/478	11mm jaw, 38mm long
SI-GI/479	14mm jaw, 48mm long
SI-GI/491	17mm jaw, 60mm long
SI-GI/492	27mm jaw, 70mm long

Plain, with hole in handle

Straight

Code	Description
SI-GI/486	11mm jaw, 38mm long
SI-GI/487	14mm jaw, 48mm long
SI-GI/488	17mm jaw, 60mm long
SI-GI/489	27mm jaw, 70mm long

Curved

Code	Description
SI-GI/483	11mm jaw, 38mm long
SI-GI/484	14mm jaw, 48mm long
SI-GI/485	17mm jaw, 60mm long
SI-GI/490	27mm jaw, 70mm long

An Overview On

Forcep Jaw Types

You will see that forceps in our catalogue will have either a traditional serrated jaw or a different design which is known as atraumatic (or non-traumatic).

Traditional Serrated Jaw



The traditional serrated jaw is generally more robust and used on forceps where the areas are not as delicate.


Atraumatic Jaw



Atraumatic jaws will be mainly found in delicate surgery areas such as vascular and cardiothoracic. The design is such, that the surface area gripped is reduced, ensuring less trauma without the loss of grip.

SPRING FORCEPS


Adsons Dissecting Forceps



127mm long

Code	Description
SI-GI/257	1-2 teeth
SI-GI/258	Lightweight, micro
SI-GI/259	Serrated
SI-GI/260	Serrated, lightweight, micro


Block End Dissecting Forceps



Serrated

Code	Description
SI-GI/255	127mm long
SI-GI/256	152mm long
SI-GI/256a	203mm long

Bonney Dissecting Forceps




178mm long

Code	Description
SI-GI/267	Serrated, 1-2 teeth
SI-GI/268	1-2 teeth

Charnley-McIndoe Dissecting Forceps



150mm long



Code	Description
SI-OR/N045a	2-3 teeth

DeBakey Atraumatic Dissecting Forceps

Straight, atraumatic jaws



1.5mm

Code	Description
SI-GI/238a	150mm long
SI-GI/240a	240mm long

2.5mm

Code	Description
SI-GI/238b	150mm long
SI-GI/240b	240mm long

2mm


Code	Description
SI-GI/238	150mm long
SI-GI/241	180mm long

2.8mm

Code	Description
SI-GI/238c	150mm long
SI-GI/240c	240mm long
SI-GI/241c	180mm long

Dissecting Forceps

Turnover end, serrated jaws



Code	Description
SI-GI/242	127mm long
SI-GI/243	152mm long

Code	Description
SI-GI/244	178mm long
SI-GI/245	203mm long

Dunhill Dissecting Forceps


Light, fine serrated jaws



Code	Description
SI-GI/457	145mm long

Fredrick Coole Dissecting Forceps

Extra fine serrated jaws




125mm long

Code	Description
SI-GI/066	Straight
SI-GI/067	Curved

SPRING FORCEPS

Gillies Dissecting Forceps


1-2 teeth



Code	Description
SI-PL/052	152mm long

McIndoe Dissecting Forceps


152mm long



Code	Description
SI-GI/476	Serrated
SI-GI/477	1-2 teeth

Jackson Burrows Dissecting Forceps


140mm long



Code	Description
SI-OR/N037	Tongue and groove, 1-2 teeth
SI-OR/N037a	Tongue and groove only

McIndoe Modified Dissecting Forceps


152mm long, 1mm tip, modified for fine powerful dissection with strong action tip



Code	Description
SI/5961	Serrated


Jean Dissecting Forceps

Serrated jaws, 140mm long



Code	Description
SI-GI/269	Fine
SI-GI/270	1-2 teeth
SI-GI/271	2-3 teeth, 178mm long

Treves Dissecting Forceps



1-2 teeth


Code	Description
SI-GI/251	127mm long
SI-GI/252	152mm long

2-3 teeth

Code	Description
SI-GI/253	127mm long
SI-GI/254	152mm long

Lanes Dissecting Forceps


1-2 teeth on serrated jaws, 140mm long



Code	Description
SI-GI/246	127mm long
SI-GI/247	152mm long
SI-GI/248	178mm long

Code	Description
SI-GI/249	2-3 teeth, 178mm long
SI-GI/250	2-3 teeth 203mm long

Waugh Dissecting Forceps



Fine, serrated

Code	Description
SI-GI/261	152mm long
SI-GI/262	178mm long
SI-GI/263	203mm long


Fine, 1-2 teeth

Code	Description
SI-GI/264	152mm long
SI-GI/265	178mm long
SI-GI/266	203mm long

TISSUE FORCEPS

Allis Tissue Forceps

Box joint



150mm long



Code	Description
SI-GI/272	3-4 teeth
SI-GI/273	4-5 teeth

190mm long

Code	Description
SI-GI/414	5-6 teeth

Duval Tissue Forceps


Box joint



Code	Description
SI-GI/279	13mm wide, 203mm long

Littlewood Tissue Forceps


Box joint, 2-3 teeth



Code	Description
SI-GI/275	190mm long

Lane Tissue Forceps


Box joint, 1-2 teeth, fenestrated



Code	Description
SI-GI/280	127mm long
SI-GI/281	152mm long
SI-GI/282	197mm long

Babcock Tissue Forceps


Box joint, serrated jaws



Code	Description
SI-GI/276	160mm long
SI-GI/277	200mm long
SI-GI/278	230mm long

Rutherford Morrison Tissue Forceps

Box joint




Code	Description
SI-GI/283	4-5 teeth, 152mm long

TISSUE FORCEPS

Shaw Tissue Forceps

Box joint




Code	Description
SI-GI/458	1-2 teeth, 178mm long



Guys Gallabin Suture-Holding Forceps


Screw joint, angled flat blades




Code	Description
SI-GI/450	125mm long
SI-GI/4402	150mm long

Russian Dissecting Forceps

With fenestrated end



Code	Description
SI-GI/SH-5503b	150mm long
SI-GI/SH-5503	200mm long
SI-GI/403	230mm long
SI-GI/SH-5503a	250mm long




Special Offer

on SSD Essentials Pack

Top up your trays with these essentials...


Halstead Mosquito Forceps

Straight/Curved
125mm (5")



Ramsey Sponge Holding Forceps

Straight, Serrated
Jaws 241mm
(9½")




Scalpel Blade Handle

No. 3




Spring Clip

Serrated Pattern,
Stainless Steel for
Fixing Slings to Leg
Splints, also for
Towel Clipping
6.3cm (2½") Wide



Spencer Wells Artery Forceps


Straight, Serrated
Jaws



SPONGE HOLDING AND STERILISING FORCEPS

Rampley Sponge Holding Forceps


Box joint, straight, serrated jaws




Code	Description
SI-GI/353	178mm long
SI-GI/354	241mm long

Lister Sinus Forceps

Box joint, serrated tip




Code	Description
SI-GI/367	127mm long
SI-GI/368	152mm long
SI-GI/369	178mm long




Foerster Sponge Holding Forceps

Box joint, straight, serrated jaws




Code	Description
SI-GI/354a	178mm long



Cheatle Sterilizer Forceps

Screw joint




Code	Description
SI-GI/364	267mm long

Harrison Bowl Steriliser Forceps

Screw joint, angled, double jaw

Code	Description
SI-GI/365	356mm long
SI-GI/366	457mm long



OTHER FORCEPS

Hunter Splinter Forceps


Serrated grooved pointed jaws



Code	Description
SI-GI/360	115mm long

Michel Forceps


Box joint, double-ended, for inserting and removing clips stainless steel



Code	Description
SI-GI/357	130mm long

Michel Skin Holding Forceps


With magazine for Michel clips



Code	Description
SI-GI/358	1-2 teeth, 127mm long

Michel Clip Extracting Forceps


Box joint



Code	Description
SI-GI/359	125mm long

Childe Approximating Forceps


With magazine for Michel clips



Code	Description
SI-GI/469	180mm long

Robins Anchoring Forceps


140mm long



Code	Description
SI-GI/453	5mm with holder for diathermy leads
SI-GI/453a	13mm with holder for suction tubing

PREMIUM QUALITY STAINLESS STEEL SCISSORS

Dressing Scissors



Blunt/Blunt

Straight

Code	Description
SI-GI/133	130mm long
SI-GI/136	146mm long
SI-GI/139	165mm long

Curved

Code	Description
SI-GI/142	130mm long
SI-GI/143	146mm long
SI-GI/144	165mm long
SI-GI/145	185mm long
SI-GI/404	147mm long

With one serrated blade

Blunt/Sharp

Straight

Code	Description
SI-GI/134	130mm long
SI-GI/137	146mm long
SI-GI/140	165mm long

Sharp/Sharp


Straight

Code	Description
SI-GI/135	130mm long
SI-GI/138	146mm long
SI-GI/141	165mm long

Curved

Code	Description
SI-GI/146	130mm long

Metzenbaum Scissors



Straight

Code	Description
SI-GI/405	180mm long
SI-GI/406	147mm long
SI-GI/406a	200mm long
SI-GI/406b	230mm long
SI-GI/406c	250mm long
SI-GI/406d	280mm long

Curved

Code	Description
SI-GI/407	185mm long
SI-GI/408	147mm long
SI-GI/408a	200mm long
SI-GI/408b	230mm long
SI-GI/408c	250mm long
SI-GI/408d	280mm long

Mayo Scissors



Straight

Code	Description
SI-GI/409	140mm long
SI-GI/409a	160mm long
SI-GI/411	190mm long
SI-GI/411a	230mm long

Curved


Code	Description
SI-GI/410	140mm long
SI-GI/410a	160mm long
SI-GI/412	190mm long
SI-GI/412a	230mm long

Dressing Scissors, Left-Handed

With one black bow, 145mm long

Blunt/Blunt

Code	Description
SI-GI/192	Straight
SI-GI/193	Curved



Blunt/Sharp

Code	Description
SI-GI/194	Straight
SI-GI/195	Curved


Sharp/Sharp

Code	Description
SI-GI/196	Straight
SI-GI/197	Curved

Littler Dissecting Scissors

Fine pattern, with eye for suture


Code	Description
SI-PL/001f	115mm long



Kilner Dissecting Scissors

Fine pattern


Code	Description
SI-PL/001a	115mm long



Fergusson Scissors


Angled on flat

Code	Description
SI-GI/148	180mm long



Paynes Suture Scissors


Code	Description
SI-GI/152	100mm long



Ligature Stitch Scissors


Fine pattern

Code	Description
SI-GI/153	115mm long



SCISSORS


Lister Bandage Scissors



Code	Description
SI-GI/154	147mm long
SI-GI/155	178mm long

Lister Bandage Scissors


Angled, with probe end



Code	Description
SI-GI/361	147mm long
SI-GI/362	178mm long

McIndoe Scissors

178mm long




Code	Description
SI-GI/149	Straight
SI-GI/150	Curved on flat
SI-GI/150sc	Supercut, curved




Thomas Strabismus Scissors

114mm long




Code	Description
SI-OP/030	Straight
SI-OP/029	Curved

Heath Stitch Scissors




Code	Description
SI-GI/151	152mm long

“Careless” Scissors



Code	Description
SI-GI/474	165mm long

Tubing Scissors



Code	Description
SI-GI/147	130mm long

Utility Scissors




Code	Description
SI-GI/466	190mm long

TUNGSTEN CARBIDE HARD INSERT SCISSORS

Dressing Scissors

Tungsten carbide hard insert, blunt/blunt



Straight


Code	Description
SI-GI/430	130mm long
SI-GI/432	145mm long
SI-GI/434	165mm long
SI-GI/436	185mm long

Curved

Code	Description
SI-GI/431	130mm long
SI-GI/433	145mm long
SI-GI/435	165mm long

Mayo Dressing Scissors

Tungsten carbide hard insert



Straight


Code	Description
SI-GI/419	140mm long
SI-GI/419a	160mm long
SI-GI/421	190mm long
SI-GI/421a	230mm long

Curved

Code	Description
SI-GI/420	140mm long
SI-GI/420a	160mm long
SI-GI/422	190mm long
SI-GI/422a	230mm long

Dressing Scissors

Tungsten carbide hard insert, blunt/sharp



Straight


Code	Description
SI-GI/423	130mm long
SI-GI/425	145mm long
SI-GI/427	165mm long
SI-GI/429	185mm long

Curved

Code	Description
SI-GI/424	130mm long
SI-GI/426	145mm long
SI-GI/428	165mm long

Mayo-Harrington Dressing Scissors

Tungsten carbide hard insert



230mm long


Code	Description
SI-GI/417	Straight
SI-GI/418	Curved

200mm long

Code	Description
SI-GI/162a	Straight
SI-GI/413b	Curved

Dressing Scissors

Tungsten carbide hard insert sharp/sharp



Straight


Code	Description
SI-GI/437	130mm long
SI-GI/439	145mm long
SI-GI/441	165mm long
SI-GI/443	185mm long

Curved

Code	Description
SI-GI/438	130mm long
SI-GI/440	145mm long
SI-GI/442	165mm long

Metzenbaum Scissors

Tungsten carbide hard insert



230mm long


Code	Description
SI-GI/445	145mm long
SI-GI/448	185mm long
SI-GI/471a	230mm long
SI-GI/471b	250mm long
SI-GI/471c	280mm long

200mm long

Code	Description
SI-GI/446	145mm long
SI-GI/449	185mm long
SI-GI/449a	200mm long
SI-GI/449b	230mm long
SI-GI/449c	250mm long
SI-GI/449d	280mm long

Baby Metzenbaum Scissors


Tungsten carbide, hard insert, 145mm long



Code	Description
SI-GI/445a	Straight
SI-GI/446a	Curved

Ligature Scissors

Tungsten carbide hard insert, curved




Code	Description
SI-GI/444	145mm long
SI-GI/447	185mm long

TUNGSTEN CARBIDE SUPERCUT SCISSORS

Mayo Harrington Scissors

Supercut, tungsten carbide



230mm long


Code	Description
SI-GI/162	Straight
SI-GI/413	Curved

200mm long

Code	Description
SI-GI/162b	Straight
SI-GI/413b	Curved

Mayo Scissors

Supercut, tungsten carbide



Straight


Code	Description
SI-GI/163	140mm long
SI-GI/163a	160mm long
SI-GI/165	190mm long
SI-GI/165a	230mm long

Curved

Code	Description
SI-GI/164	140mm long
SI-GI/164a	160mm long
SI-GI/166	190mm long
SI-GI/166a	230mm long

Metzenbaum Dissecting Scissors

Supercut, tungsten carbide



Straight


Code	Description
SI-GI/188tc	185mm long
SI-GI/188tca	200mm long
SI-GI/188tcb	230mm long
SI-GI/188tcc	250mm long
SI-GI/188tcd	280mm long

Curved

Code	Description
SI-GI/191tc	147mm
SI-GI/189tc	185mm
SI-GI/186tca	200mm
SI-GI/188tcb	230mm
SI-GI/188tcc	250mm
SI-GI/188tcd	280mm

Dressing Scissors

Supercut, tungsten carbide



Blunt/Blunt

Straight

Code	Description
SI-GI/174	130mm long
SI-GI/176	145mm long
SI-GI/178	165mm long
SI-GI/180	185mm long

Curved

Code	Description
SI-GI/175	130mm long
SI-GI/177	145mm long
SI-GI/179	165mm long
SI-GI/180a	185mm long

Blunt/Sharp

Straight

Code	Description
SI-GI/167	130mm long
SI-GI/169	145mm long
SI-GI/171	165mm long
SI-GI/173	185mm long

Curved

Code	Description
SI-GI/168	130mm long
SI-GI/170	145mm long
SI-GI/172	165mm long
SI-GI/173a	185mm long

Sharp/Sharp

Straight

Code	Description
SI-GI/181	130mm long
SI-GI/183	145mm long
SI-GI/185	165mm long
SI-GI/187	185mm long

Curved

Code	Description
SI-GI/182	130mm long
SI-GI/184	145mm long
SI-GI/186	165mm long
SI-GI/187a	185mm long

Guide for Gold and Black Bows

It's important to know the difference for the intended use of each scissor:



2 Gold Bows
Tungsten carbide insert or padding for a scissor



2 Black Bows
Supercut edge



1 Black and 1 Gold Bow
Supercut tungsten carbide




1 Black Bow (left ring)
Left-handed instrument

SUPERCUT OPERATING SCISSORS

Dressing Scissors

Supercut, blunt/blunt




145mm long

Code	Description
SI-GI/160	Straight
SI-GI/161	Curved

Dressing Scissors

Supercut, blunt/sharp



145mm long

Code	Description
SI-GI/156	Straight
SI-GI/157	Curved

Dressing Scissors

Supercut, sharp/sharp




145mm long

Code	Description
SI-GI/158	Straight
SI-GI/159	Curved

Mayo Scissors

Supercut



Straight


Code	Description
SI-GI/166a	140mm long
SI-GI/166e	160mm long
SI-GI/166b	170mm long

Curved

Code	Description
SI-GI/166c	140mm long
SI-GI/166f	160mm long
SI-GI/166d	170mm long

Boyd Scissors


Supercut, blunt/blunt



Code	Description
SI-THO/047	180mm long

Metzenbaum Dissecting Scissors

Supercut



Straight (length)


Code	Description
SI-GI/190	147mm
SI-GI/188	185mm
SI-GI/188a	200mm
SI-GI/188b	230mm
SI-GI/188c	250mm
SI-GI/188d	280mm

Curved (length)

Code	Description
SI-GI/191	147mm
SI-GI/189	185mm
SI-GI/186a	200mm
SI-GI/186b	230mm
SI-GI/186c	250mm
SI-GI/186d	280mm

Stevens Tenotomy Scissors

Supercut



Curved/blunt

Code	Description
SI-OP/424	114mm long
SI-OP/424a	125mm long
SI-OP/424gb	135mm long

Curved/sharp

Code	Description
SI-OP/425	114mm long
SI-OP/425a	140mm long
SI-OP/425b	160mm long


With cutout, 110mm long

Code	Description
SI-OP/426	Straight/blunt, 114mm long
SI-OP/427	Straight/sharp, 114mm long
SI-OP/427c	Straight/sharp, 180mm long

Code	Description
SI-OP/427a	Straight/blunt
SI-OP/427b	Curved/blunt

Strabismus Scissors

Supercut




114mm long

Code	Description
SI-GI/475	Straight
SI-GI/475a	Curved

Iris Scissors

Supercut




114mm long

Code	Description
SI-OP/476	Straight
SI-OP/476a	Curved

NEEDLE HOLDERS

Bruce Clarke Needle Holder


Screw joint



Code	Description
SI-GI/205	127mm long

Lawrence Needle Holder


Fine jaw, 150mm long



Code	Description
SI-PL/012a	
SI-PL/012b	Tungsten carbide jaw

Crile-Wood Needle Holder


Straight, baby, fine



Code	Description
SI-GI/199	150mm long

Needle Holder


Surgical Holdings pattern, box joint, lightweight, with curved plain jaws and one offset bow for thumb



Code	Description
SI-GI/209	152mm long

Kilner Needle Holder


Box joint, cross-serrated jaw



Code	Description
SI-GI/200	165mm long
SI-GI/201	133mm long

MacPhail Needle Holder

Screw joint, with one copper jaw




Code	Description
SI-GI/206	159mm long
SI-GI/207	178mm long
SI-GI/208	209mm long

NEEDLE HOLDERS

Mathieu Needle Holder


Screw joint



Code	Description
SI-GI/210	172mm long
SI-GI/211	209mm long
SI-GI/212	243mm long


Naughton Morgan Needle Holder

Box joints, cross serrated jaws



Mayo Needle Holder

Box joint, straight



Wide cross serrated jaws

Code	Description
SI-GI/224	159mm long
SI-GI/398	203mm long

Narrow cross serrated jaws

Code	Description
SI-GI/225	159mm long


Tungsten carbide

Code	Description
SI-GI/399	185mm long
SI-GI/400	254mm long

Code	Description
SI-GI/401	185mm long
SI-GI/402	254mm long

Rochester-Mayo Needle Holder

Screw joint



Code	Description
SI-GI/213	178mm long

Choosing the most suitable needle size

Our tungsten carbide inserts are available in a variety of styles. Choosing the most suitable needle size for each procedure helps give the surgeon maximum grip and control.

Coarse 0.65mm

For powerful needle holder patterns and wire twisters. Maximum holding capabilities for heavy needles or wire.



Standard 0.5mm

Used, for example, in Mayo-Hegar, Mathieu Needle Holders. For needles 3-0 and larger.



Fine 0.4mm

For delicate needle holders such as Crile-Wood and DeBakey. For needles and sutures between 4-0 and 6-0.



Micro 0.2-0.3mm

For micro needle holders such as Castroviejo, Microvascular and Mini-Ryder. For delicate needles and sutures between 7-0 and 10-0.



Smooth


Smooth tungsten carbide surface reduces damage to delicate needles and suture of 6-0 and finer.



TUNGSTEN CARBIDE NEEDLE HOLDERS

Berry Needle Holder/Wire Twister


Tungsten carbide jaws



Code	Description
SI-GI/467	200mm long

Finochietto Needle Holder


Tungsten carbide jaws



Code	Description
SI-GI/222	215mm long
SI-GI/223	265mm long

Crile-Wood Needle Holder


Tungsten carbide insert, box joint, satin finish, straight



Code	Description
SI-GI/217	159mm long
SI-GI/218	178mm long
SI-GI/198	150mm long

Halsey Needle Holder

Tungsten carbide insert, box joint, satin finish, straight



Code	Description
SI-GI/220	127mm long

Derf Needle Holder


Tungsten carbide insert



Code	Description
SI-GI/BUPHF002	115mm long

Kilner Needle Holder

Tungsten carbide insert




Code	Description
SI-GI/221	140mm long

TUNGSTEN CARBIDE NEEDLE HOLDERS

Masing Needle Holder


Tungsten carbide insert, one angled bow



Code	Description
SI-GI/009R	130mm long

Nievert Needle Holder


Tungsten carbide insert, box joint, satin finish, straight



Code	Description
SI-GI/219	127mm long

DeBakey Needle Holder


250mm long



Code
THO/ SI-5642a

Ryder Needle Holder

Tungsten carbide insert, box joint, satin finish, straight



0.2mm jaw pitch


Code	Description
SI-GI/SH-5197a	135mm long
SI-GI/SH-5197b	155mm long
SI-GI/SH-5197c	175mm long
SI-GI/SH-5197d	210mm long

0.4mm jaw pitch

Code	Description
SI-GI/SH-5198a	135mm long
SI-GI/SH-5198b	155mm long
SI-GI/SH-5198c	175mm long
SI-GI/SH-5198d	210mm long

Sarot Needle Holder

Diamond jaw, tungsten carbide, gold plated bows



Code	Description
SI-GI/451	181mm long

Tungsten Carbide Needle Holders

Features and Benefits

Jaw

Tungsten carbide inserts bonded into the jaws for a long lasting, secure grip

Bevelled Edges

Prevents suture from snagging

Laser Etching

Date and catalogue number are etched for easy re-ordering. Laser etch will not harm instrument like traditional methods

Satin Finish

Reduces glare

Ratchets

Precisely cut for effortless locking/unlocking

Teeth

Fine, hard teeth, diamond cut into an insert of hardened tungsten carbide for superior needle grasping and suture tying capabilities

Vacuum Brazing

Combines bonding of the tungsten carbide insert with the simultaneous vacuum hardening of the instrument for greater durability and longevity

Shanks


Precision balanced for comfort and control

Gold Handles

For easy identification

Mayo-Hegar Needle Holder

Tungsten carbide insert, box joint, satin finish, straight




Code	Description
SI-GI/472	130mm long
SI-GI/214	159mm long
SI-GI/215	178mm long
SI-GI/216	209mm long

LONGER SIZES AVAILABLE

COMBINED SCISSOR AND NEEDLE HOLDERS

Olsen-Hegar Scissor/Needle Holder


With tungsten carbide insert



Code	Description
SI-GI/229	147mm long
SI-GI/230	165mm long

Olsen-Hegar Scissor/Needle Holder


Screw joint



Code	Description
SI-GI/226	140mm long
SI-GI/227	165mm long
SI-GI/228	190mm long

Gillies Scissor/Needle Holder

Right-hand, with tungsten carbide insert



Code	Description
SI-GI/235	158mm long

Gillies Scissor/Needle Holder


Heavy pattern, 158mm long



Code	Description
SI-GI/233	Right-hand
SI-GI/234	Left-hand

Gillies Scissor/Needle Holder


158mm long



Code	Description
SI-GI/231	Right-hand
SI-GI/232	Left-hand

Foster Gillies Scissor/Needle Holder

Fine pattern




Code	Description
SI-GI/236	127mm long

RONGEURS

Beck Rongeur


Straight, for intervertebral disc



Code	Description
SI-NEU/125	2mm bite, 130mm shaft

Pennybacker Rongeur

120mm shaft




Code	Description
SI-NEU/133	3mm bite
SI-NEU/107	4mm bite

Cushing Rongeur

130mm shaft




Code	Description
SI-NEU/127-2	2mm bite
SI-NEU/127-3	3mm bite
SI-NEU/127-4	4mm bite
SI-NEU/127-5	5mm bite
SI-NEU/127-6	6mm bite



St Bart's Rongeur


114mm shaft



Code	Description
SI-ENT/140a	3mm bite

Cushing Rongeur

130mm shaft



150° upward cutting


Code	Description
SI-NEU/128c	2mm bite
SI-NEU/128	3mm bite
SI-NEU/128a	4mm bite
SI-NEU/128b	5mm bite

150° downward cutting

Code	Description
SI-NEU/130c	2mm bite
SI-NEU/130	3mm bite
SI-NEU/130a	4mm bite
SI-NEU/130b	5mm bite

St Bart's Rongeur

140°, 3mm, 114mm shaft




Code	Description
SI-ENT/140b	Upward cutting
SI-ENT/140ac	Downward cutting

BAG AND TOWEL CLIPS

Spring Clip


Serrated pattern, stainless steel, for fixing slings to leg splints, also for towel clipping



Code	Description
SI-GI/344	63mm wide

Spring Bag Clip


Serrated



Code	Description
SI-GI/345	26mm wide

Jones Towel Clip


Superior quality spring forcep type



Code	Description
SI-GI/346	89mm long

Shardles Towel Clip


Spring curved, cross action



Code	Description
SI-GI/347	77mm long

Backhaus Towel Clip


Superior quality, box joint



Code	Description
SI-GI/348	89mm long
SI-GI/348a	102mm long
SI-GI/349	130mm long

Moynihan Towel Clip


Screw joint, 2-2 teeth, for sides of incision



Code	Description
SI-GI/470	191mm long

Ball and Socket Towel Forceps

Box joint




Code	Description
SI-GI/350	89mm long
SI-GI/351	115mm long
SI-GI/352	140mm long

RETRACTORS

Goligher Retractor

Complete with 2 centre blades:
1 x 90mm wide x 50mm deep
1 x 90mm wide x 90mm deep

and 2 pairs of lateral blades:
1 pair 25mm wide x 45mm deep
1 pair 50mm wide x 90mm deep



Code

SI-GI/305

Balfour Abdominal Retractor

Trivalve, frame 203mm

Child


102mm opening, 51mm deep lateral blades,
38mm deep centre blade

Adult

178mm opening, 64mm deep lateral blades,
44mm deep centre blade

Extra large

178mm opening, 102mm deep lateral blades,
108mm deep centre blade




Code	Description
SI-GI/373	Child
SI-GI/374	Adult
SI-GI/375	Extra large

Miles Retractor

Complete with 2 centre blades:
1 x 90mm wide x 50mm deep
1 x 90mm wide x 75mm deep

and 2 pairs of lateral blades:
1 pair 80mm wide x 35mm deep
1 pair 80mm wide x 75mm deep




Code	Description
SI-GI/376	340mm long

OTHER BLADE SIZES AVAILABLE

Balfour Abdominal Retractor


Bivalve, baby



Code	Description
SI-GI/372	Fully open 89mm

Berkeley Abdominal Retractor

Parallel action




Code	Description
SI-GI/461	Blades 35mm x 60mm deep, 215mm long

TURNER WARWICK RETRACTORS

Turner Warwick Universal Self Retaining Retractor

Code	Description
SI-GI/520	With frames A/B/C/D and blades A/B/C/D/E/F/K/L
SI-GI/521-A	Circular frame A only 178mm
SI-GI/521-B	Circular frame B only 229mm
SI-GI/521-C	Circular frame C only 305mm
SI-GI/521-D	Oval frame D only 475mm x 330mm
SI-GI/521-E	Oval frame E only 370mm x 300mm



Turner Warwick Needleholder with Tungsten Carbide Jaws

Curved shanks with one short bow


Code	Description
SI-5835	230mm long
SI-5835A	200mm long

Blade A




Code	Description
SI-GI/522-A	Blade only - 19mm x 76mm

Blade B




Code	Description
SI-GI/522-B	Blade only - 38mm x 102mm

Blade C



Code	Description
SI-GI/522-C	Blade only - 57mm x 89mm


Blade D



Code	Description
SI-GI/522-D	Blade only - 57mm x 114mm


TURNER WARWICK RETRACTORS

Blade E




Code	Description
SI-GI/522-E	Blade only - curved stainless steel blade 51mm x 38mm

Blade F




Code	Description
SI-GI/522-F	Blade only - curved stainless steel blade 76mm x 60mm

Blade G




Code	Description
SI-GI/522-G	Blade only - flat stainless steel blade 89mm x 60mm

Blade J




Code	Description
SI-GI/522-J	Blade only - single hook, 13mm deep

Blade K



Code	Description
SI SC-GI/522-K	Blade only - flat with clip 55mm x 90mm

Blade L




Code	Description
SI-GI/522-L	Blade only - flat with clip 55mm x 120mm

RETRACTORS

Collin Retractor


110mm fully open, 235mm long



Code	Description
SI-GI/322	Complete with centre blade 55 x 40mm and one pair of blades 38mm x 60mm
SI-GI/322a	Without centre blade

Gelpi Retractor


Self retaining



Code	Description
SI-GY/013a	140mm long
SI-GY/013b	178mm long


Joll Thyroid Retractor

Self retaining, with lock-in nut




Code	Description
SI-GI/363	152mm long

Gossett Abdominal Retractor




Code	Description
SI-GI/303	Large, fully open 185mm long
SI-GI/304	Small, fully open 95mm long

Pozzi Abdominal Retractors



Code	Description
SI-GI/370	Swivel blades: 38 x 60mm, bivalve, 254mm long




Code	Description
SI-GI/371	Swivel blades: 38 x 59mm, trivalve, 203mm long

RETRACTORS

Bowlby Dilator


Curved



Code	Description
SI-GI/495	145mm long

Travers Retractor


Self retaining with internal rack, 210mm long



Code	Description
SI-GI/323	4 x 5 teeth, 95mm fully open

Irwin Retractor


Self retaining, with cam rack, 4 x 4 teeth



Code	Description
SI-GI/324	178mm long

West Retractor


Self retaining, 140mm long



Code	Description
SI-GI/326	3 x 4 teeth, 60mm fully open
SI-GI/326cr	3 x 4 teeth, 60mm fully open, with cam rack

Norfolk and Norwich Retractor


4 x 5 teeth 50mm blades x 223mm long



Code	Description
SI-GI/325	Standard rack

Weitlander Retractor

3-4 teeth



135mm long


Code	Description
SI-GI/496	Blunt
SI-GI/496a	Sharp
SI-GI/496b	Sharp, 2-3 teeth
SI-GI/496c	Blunt, 2-3 teeth

165mm long

Code	Description
SI-GI/497	Blunt
SI-GI/497a	Sharp

RETRACTORS

Kocher Retractor



230mm long


Code	Description
SI-GI/378	36 x 20mm blade

250mm long

Code	Description
SI-GI/378a	67 x 25mm blade
SI-GI/378b	80 x 42mm blade
SI-GI/378c	80 x 50mm blade
SI-GI/378d	80 x 60mm blade

Kocher Langenbeck Retractor

215mm long




Code	Description
SI-GI/336sp	25 x 6mm blade
SI-GI/379a	35 x 8mm blade
SI-GI/379b	35 x 11mm blade
SI-GI/379c	35 x 15mm blade
SI-GI/357sp	41 x 11mm blade

Code	Description
SI-THO/cutsp015	55 x 11mm blade
SI-GI/379d	70 x 15mm blade
SI-GI/379e	80 x 12mm blade
SI-GI/379f	80 x 16mm blade


Langenbeck Retractor

Single ended



Code	Description
SI-GI/336	Blade 22mm x 6mm
SI-GI/337	Blade 45mm x 13mm
SI-GI/337a	Blade 65mm x 13mm
SI-GI/338a	Blade 75mm x 25mm
SI-GI/338a1	Blade 65mm x 25mm

Landon Retractor




Code	Description
SI-GI/494	230mm



Canny Ryall Navy Retractor

Double-ended, with one single hook blade 5mm wide, and one fenestrated blade



Code	Description
SI-GI/343	Fenestrated blade 16mm wide, 190mm long
SI-GI/343a	Fenestrated blade 19mm wide, 190mm long
SI-GI/343b	Fenestrated blade 25mm wide, 220mm long
SI-GI/343c	Fenestrated blade 30mm wide, 220mm long

Canny Ryall Army Retractor

Double prong, 13mm wide, 165mm long




Code	Description
SI-GI/342	Large, 19mm wide
SI-GI/342a	Small, 11mm wide



Czerny Retractor


Double ended



Code	Description
SI-GI/341	171mm long

Lane Retractor

Double ended blade, 229mm long




Code	Description
SI-GI/339	13mm x 19mm wide
SI-GI/340	25mm x 32mm wide

RETRACTORS

Morris Retractor


Single blade, 229mm long



Code	Description
SI-GI/330	38 x 51mm
SI-GI/331	51 x 51mm
SI-GI/332	64 x 51mm


Morris Retractor

Double ended, 216mm long




Code	Description
SI-GI/327	32 x 51mm and 38 x 51mm
SI-GI/328	44 x 51mm and 51 x 51mm
SI-GI/329	64 x 51mm and 38 x 51mm

Moyniham Retractor




Code	Description
SI-GI/462	Single ended, blade 13 x 22mm, 215mm long

Deaver Abdominal Retractors



Code	Description
SI-GI/306	Set of 5, 315mm overall length
SI-GI/307	25mm wide
SI-GI/308	25mm wide
SI-GI/309	38mm wide
SI-GI/310	51mm wide
SI-GI/311	76mm wide


Elliot (Modified Deaver) Retractor



Code	Description
SI-GI/308M	With flat handle

Dyball and Lock Retractor

For liver




Code	Description
SI-GI/335	64 x 144mm, 241mm long

OTHER SIZES AVAILABLE - PRICE ON APPLICATION

Deaver Abdominal Retractors

With round handles



Code	Description
SI-GI/312	Set of 5, 305mm overall length
SI-GI/313	25mm wide
SI-GI/314	25mm wide
SI-GI/315	38mm wide
SI-GI/316	51mm wide
SI-GI/317	76mm wide

Kelly Retractor Blade

267mm long




Code	Description
SI-GI/333	38 x 191mm
SI-GI/334	57 x 191mm

RETRACTORS

Durham Retractor


With fenestrated blade, 216mm long



Code	Description
SI-GI/300	10mm wide
SI-GI/301	14mm wide
SI-GI/302	20mm wide

Ballance (Volkman) Retractor

203mm long



Sharp

Code	Description
SI-GI/463	1 prong
SI-GI/463a	2 prong
SI-GI/463b	3 prong
SI-GI/463c	4 prong
SI-GI/463d	6 prong

Blunt

Code	Description
SI-GI/463e	1 prong
SI-GI/463f	2 prong
SI-GI/463g	3 prong
SI-GI/463h	4 prong
SI-GI/463i	6 prong

Volkman Retractor



Code	Description
SI-OR/400	Sharp, 220mm long
SI-OR/400a	Blunt, 220mm long
SI-OR/400b	Sharp, 245mm long

Kilner Retractor

Double ended, cat's paw



Code	Description
SI-PL/056	152mm long

Kocher Retractor



Code	Description
SI-OR/N036c	Sharp, 205mm long
SI-OR/N036	Sharp, 220mm long
SI-OR/N036b	Blunt, 220mm long

Senn-Miller Retractor

170mm long



Code
SI-SH/5881

Single Hook Retractor

Sharp



Code	Description
SI-GI/SH5671	152mm long
SI-PL/065b	Half size, fine

Gillies Skin Hook

Solid handle, 190mm long



Code	Description
SI-PL/065a	Standard, fine
SI-PL/065b	Half size, fine

DISSECTORS AND PROBES

MacDonald Dissector Double Ended



Code	Description
SI-GI/285	191mm long

Brodie Director Probe Ended, Silver Plated



Code	Description
SI-GI/286	165mm long

Watson Cheyne Dissector

Fine, double-ended, with probe end



Code	Description
SI-GI/287	127mm long
SI-GI/288	191mm long
SI-GI/288a	230mm long
SI-GI/288b	255mm long

Probe, Malleable Silver Plated



Code	Description
SI-GI/289a	127mm long
SI-GI/289	152mm long
SI-GI/290	178mm long
SI-GI/291	203mm long
SI-GI/292	255mm long
SI-GI/292a	305mm long

Probe, Sinus Silver



Code	Description
SI-GI/464	152mm long

Probe, Fish Tail Silver Plated



Code	Description
SI-GI/293	127mm long
SI-GI/294	152mm long
SI-GI/295	203mm long
SI-GI/296	254mm long

Jobson Horne Probe Dressing (Aural)


With one flat tapered end, serrated edge and one small ring end, stainless steel



Code	Description
SI-GI/297	140mm long
SI-GI/298	178mm long

Bookwalter® Retractor Kit III

A table-fixed retractor system designed by surgeons to efficiently and reliably hold surgical instruments in place during surgery while allowing greater visibility.



Kit includes Bookwalter® Kit:

Table Post and Coupling, Horizontal Flex Bar, Segmented Rings, Tilt Ratchet Mechanism, Kelly Retractors, Blades and Ratchets

MISCELLANEOUS


Aneurysm Needle

Small, for infusions




Code	Description
SI-GI/472a	146mm long

Syme Aneurysm Needle




Code	Description
SI-GI/473	165mm long

Ring Cutter



Code	Description
SI-GI/452	Circular saw blade, 165mm long
SI-GI/452a	Spare blade

Mayo Instrument Pin



Code	Description
SI-GI/355	102mm long
SI-GI/356	152mm long

Baron Scalpel Handle


130mm long



Code	Description
SI-PM/003	No 3
SI-PM/004	No 4

Scalpel Handle

127mm long



Code	Description
SI-PM/007	No 3
SI-PM/008	No 4
SI-PM/008a	No 4 L
SI-PM/001a	No 7

Endoscope Re-fibre

Surgical Holdings manufacture their own fibre bundles in house. We cover small and larger diameter scopes and only use high quality German fibre optic material and stainless steel tubing.

Quality components, globally sourced:

- ✓ Medical grade stainless steel tubing
- ✓ German engineered new glass fibre
- ✓ New fibre cone as standard
- ✓ New solderable sapphire distal window
- ✓ New light post colour ring
- ✓ New objective in 95% of cases
- ✓ Replacement 'O' ring to prevent leakage



After repair




Before repair

HOLLOW WARE

Kidney Dish


Stainless steel



Code	Description
SI-GI/380	152mm
SI-GI/381	203mm
SI-GI/382	254mm
SI-GI/383	305mm

Lotion Bowl


Without lid, stainless steel



Code	Description
SI-GI/384	102mm
SI-GI/385	152mm
SI-GI/386	203mm
SI-GI/387	254mm
SI-GI/388	305mm


Instrument Tray

Stainless steel



Code	Description
SI-GI/389	254 x 203 x 50mm
SI-GI/390	Lid only for 254 x 203 x 50mm
SI-GI/391	305 x 254 x 50mm
SI-GI/392	Lid only for 305 x 254 x 50mm
SI-GI/393	356 x 305 x 50mm
SI-GI/394	Lid only for 356 x 305 x 50mm

Cheate Jar



Code	Description
SI-GI/397	Stainless steel, 203mm x 50mm

Measuring Jug




Code	Description
SI-GI/395	1 litre capacity

Petri Dish



Code	Description
SI-SH/002149	100mm, stainless steel, with lid

Gallipot




Code	Description
SI-GI/396	Stainless steel, 63mm

SILICONE MATTING & STERILISING BASKETS

Silicone Matting

Medical grade for storage, sterilisation and transportation of instruments.


Cushion, protect, and organise delicate surgical instruments during storage, transportation, and sterilisation. Perforations allow for ventilation during the sterilisation process. The silicone mat can be cut to fit any size container.



502mm long x 242mm wide
7mm holes, 10mm between holes

Code


SI/1004M



515mm long x 240mm wide
3mm holes, 6mm between holes

Code

SI/1001M



570mm long x 550mm wide
4mm holes, 13mm between holes

Code

SI/1002M

Sterilising Baskets

DIN Basket

Wire mesh base with perforated sides and two carrying handles

Code	Description
SI-GI/CO006	230 x 130 x 60mm
SI-GI/CO007	240 x 130 x 50mm
SI-GI/CO008	240 x 130 x 60mm
SI-GI/CO009	240 x 250 x 50mm
SI-GI/CO010	240 x 250 x 60mm
SI-GI/CO014	300 x 300 x 50mm
SI-GI/CO015	350 x 350 x 60mm
SI-GI/CO011	480 x 250 x 50mm
SI-GI/CO012	480 x 250 x 60mm

Code	Description
SI-GI/CO013	480 x 250 x 70mm
SI-GI/CO016	480 x 400 x 60mm
SI-GI/CO017	480 x 480 x 50mm
SI-GI/CO018	480 x 480 x 60mm
SI-GI/CO019	500 x 300 x 50mm
SI-GI/CO020	550 x 430 x 50mm
SI-GI/CO021	600 x 300 x 60mm



DIN Basket

Wire mesh basket with two carrying handles

Code	Description
SI-GI/CO022	240 x 250 x 50mm
SI-GI/CO023	240 x 250 x 60mm
SI-GI/CO024	240 x 250 x 70mm
SI-GI/CO031	300 x 300 x 60mm
SI-GI/CO025	480 x 250 x 50mm

Code	Description
SI-GI/CO026	480 x 250 x 60mm
SI-GI/CO027	480 x 250 x 70mm
SI-GI/CO028	480 x 250 x 100mm
SI-GI/CO029	480 x 480 x 50mm
SI-GI/CO030	480 x 480 x 70mm

DIN Basket Lids

For perforated and wire mesh baskets

Code	Description
SI-GI/CO032	240 x 250mm
SI-GI/CO033	300 x 300mm
SI-GI/CO034	480 x 250mm
SI-GI/CO035	480 x 480mm

Did you know...


Organising and fixing instruments in place can aid the cleaning process.

Help keep your instruments in good condition for longer with our baskets and trays. Available in custom sizes, to suit your requirements perfectly.

STERILISING BASKETS

ENT Basket

Wire mesh base with perforated sides and two carrying handles, supplied with lid



Code	Description
SI-GI/CO036	480 x 250 x 160mm
SI-GI/CO037	Fixing clamps

Arthroscope Basket

Wire mesh base with two fixing clamps



Code	Description
SI-GI/CO001	290 x 80 x 52mm

Endoscope Basket


With small lens basket

Code	Description
SI-GI/CO002	6400 x 150 x 77mm

Wire mesh base with two fixing clamps


Code	Description
SI-GI/CO003	S - 460 x 80 x 52mm
SI-GI/CO004	D - 480 x 80 x 52mm
SI-GI/CO005	D - 680 x 80 x 52mm

Sieve



Code	Description
SI-GI/5071	203mm x 50mm ø

Symmetry Surgical Quadlock Containers



Surgical Holdings offer the unique Quadlock container. The Patented design offers 4 locks, as opposed to the traditional 2, to ensure extra security and peace of mind.

Fine Mesh Basket



With hinged take off / lockable lid

Code	Description
SI-GI/CO041	120 x 155 x 40mm
SI-GI/CO044	120 x 155 x 50mm
SI-GI/CO042	235 x 155 x 40mm
SI-GI/CO045	235 x 155 x 50mm
SI-GI/CO043	375 x 260 x 40mm

With separate lid

Code	Description
SI-GI/CO038	120 x 155 x 40mm
SI-GI/CO039	235 x 155 x 40mm
SI-GI/CO040	375 x 260 x 40mm

OPTICAL FIBRE CABLES

Our optical fibre cables are compatible with all major scopes and illuminators. Every light lead is high temperature xenon compatible, rated to the highest level of thermal resistance and tested to ensure the highest light output possible.

Choice of 10 colours for easy identification

Choice of lengths and sheathing types



Free on-site light lead testing

We'll come and go through your lead stock to ensure your leads are fit for use



Reinforced for longer life

How to order

Select an option from each of the following steps to create your part number.

Part number =


Step 1 Step 2 Step 3 Step 4 Step 5 Step 6 Step 7

Example: 474810180FRGY is a Storz Proximal, Storz Distal, 4.8mm Bundle, 254mm Length, Fully Reinforced Sheathing in Grey

Step 1: Proximal

All end tips are stainless steel


ACMI



Ø 0.299

1

Wolf/Dyonics



Ø 0.354

2


Long Wolf



Ø 0.354

L2

Olympus/ACMI



Retaining ring


3

Storz



4

Pilling



5

Universal



3/8-24 or 5/16-24 thread

U

OPTICAL FIBRE CABLES


Step 2: Distal

All end tips are stainless steel

Wolf Snap On


6

Ø 0.355



Storz


7



ACMI Snap-On


8

Ø 0.302



Olympus


9



Wolf Clip-On


10

Ø 0.355



Pilling


11



ACMI Clip-On


12

Ø 0.302



Mini Olympus


13



2.8 mm & 4.25 mm only

Universal

U




3/8-24 or 5/16-24 thread

Step 3: Bundle Size

Place gauge on bundle to determine diameter


2.2mm

22



2.8mm

28



Mini Olympus only


3.0mm

30




3.8mm

38



4.25mm


425



Mini Olympus only

4.8mm

48



Step 4: Length

6ft

6

8ft

8

10ft

10

Call us for custom lengths

Step 5: Connector Angle

180°

180

45°

45

90°

90

Step 6: Sheathing Type

Non-reinforced

0

6" Reinforced*

6

Fully Reinforced*

FR

*Reinforced options are not available with clear sheathing

Step 7: Sheathing Colour

Black



BK

Blue



BL

Clear



CL

Green



GR

Gray




GY

Red



RD

Pink



PK

White




WH

Yellow



YL

4.8mm only



Contact Information



www.shark-instruments.com



info@shark-instruments.com



052-4598813



Mohallah Chah Tarkhana Nargate Street
Jinnah Park, Colony no.2, Sialkot.