




2025 INTESTINAL CATALOG

www.shark-instruments.com

Intestinal


Forceps	01
Clamps	03
Miscellaneous	04



FORCEPS

Desjardin Gall Stone Forceps


Screw joint



Code	Description
SI-IN/001a	Small, 220mm long
SI-IN/001b	Large, 238mm long

Lane Stomach Holding Forceps


Straight, serrated jaws, screw joint



Code	Description
SI-IN/008	203mm long

Babcock Intestinal Forceps


Straight, extra light



Code	Description
SI-GI/278	203mm long

Lahey Cholecystectomy Forceps


Curved, longitudinal serrations, heavier pattern



Code	Description
SI-IN/011	76mm blades, 197mm long
SI-IN/011a	225mm long


Maingot Intestinal Forceps

Straight, serrated tip, fenestrated grip




Code	Description
SI-IN/006	203mm long

Lahey Sweet Forceps



Code	Description
SI-IN/011b	200mm long
SI-IN/011c	225mm long
SI-IN/011d	230mm long



Denis Browne Intestinal Forceps

Straight, serrated tips, fenestrated jaws




Code	Description
SI-IN/007	178mm long

FORCEPS

Moynihan Cholecystectomy Forceps


Curved on flat, screw joint



Code	Description
SI-IN/012	215mm long

Lane Intestinal Forceps


110mm blades, longitudinal serrations
240mm long




Code	Description
SI-IN/013	Straight
SI-IN/013b	Curved

Kocher Intestinal Forceps

With light spring blades and longitudinal serrations




Code	Description
SI-IN/015a	Straight, 90mm blades, 210mm long
SI-IN/015b	Curved, 85mm blades, 210mm long



Code	Description
SI-IN/015c	Straight, 120mm blades, 250mm long
SI-IN/015d	Curved, 1150mm blades, 250mm long
SI-IN/015e	Straight, 150mm blades, 250mm long
SI-IN/015f	Curved, 145mm blades, 250mm long

Seton Pringle Intestinal Forceps


Straight, longitudinal serrations
178mm long




Code	Description
SI-IN/014	Fine blades 2.5mm wide

FORCEPS AND CLAMPS

Lang Stevenson Intestinal Forceps



Code	Description
SI-IN/018	82mm jaw, keyed blades, 203mm long



Parker Kerr Intestinal Forceps/Clamp

Blades 89mm, 254mm long


Code	Description
SI-IN/017a	Curved
SI-IN/017c	Straight



Guard/Shoe


SC-IN/017b curved, fenestrated inner or outer curve, 89mm

Code	Description
SI-IN/017d	Straight, 89mm



Schoemaker-Ogilvie Colectomy Forceps


Straight, hollow blades




Code	Description
SI-IN/020	171mm long
SI-IN/020a	203mm long

Hartman Intestinal Clamp

With diagonal serrations, 200mm long




Code	Description
SI-IN/015g	Straight
SI-IN/015h	Curved



Payr Intestinal Clamp

Crushing, lever action, longitudinal serrations




Code	Description
SI-IN/003	Extra small, 50mm blades, 150mm long
SI-IN/002a	Small, 65mm blades, 210mm long
SI-IN/002b	Medium, 100mm blades, 240mm long
SI-IN/002c	Large, 150mm blades graduated in inches, 360mm long

CLAMPS AND MISCELLANEOUS

Doyen Intestinal Clamp


240mm long



Code	Description
SI-IN/004a	Straight
SI-IN/004b	Curved


Lane Twin Stomach Clamp

Longitudinal serrations
large, 150mm blades, 305mm long



Code	Description
SI-IN/009b	Straight
SI-IN/009a	Curved

Zachary Cope Anastomosis Clamp



Code	Description
SI-IN/010a	Set: 3 clamps, lever action closing device and tommy bar, 233mm long
SI-IN/010b	Set of 3 clamps
SI-IN/010c	Tommy bar for set, 103mm long

Bakes Bile Duct Dilators

Malleable stems, 203mm long



Taking care of your instruments

It is really important to ensure your instruments remain safe.

To do this, you will need a device which allows you to magnify distal areas of the instruments to check for:


- Burrs
- Corrosion
- Cracks




Code	Description
SI-IN/016a	Set of nine
SI-IN/016c	3mm ø head
SI-IN/016d	4mm ø head
SI-IN/016e	5mm ø head
SI-IN/016f	6mm ø head

Code	Description
SI-IN/016g	7mm ø head
SI-IN/016h	8mm ø head
SI-IN/016i	9mm ø head
SI-IN/016j	10mm ø head
SI-IN/016k	11mm ø head

Yankauer Suction Tube



Code	Description
SI-IN/019	305mm long



Contact Information



www.shark-instruments.com



info@shark-instruments.com



052-4598813



Mohallah Chah Tarkhana Nargate Street
Jinnah Park, Colony no.2, Sialkot.